


PhO₂x Box

Culture as nature
intended



 Baker



PhO₂x Box gas-controlled cell culture chamber

Entry level solution for physoxia and hypoxia experiments

The PhO₂x Box system is comprised of a Gas Controller unit (offering O₂ and CO₂ control) and a Cell Culture Chamber.

The Cell Culture Chamber is placed inside any incubator or workstation for temperature control, or on a lab bench for room temperature applications. PhO₂x Box offers experimental flexibility and precision atmosphere within your existing laboratory infrastructure.

Baker's PhO₂x Box is a novel, easy to use and economical Physoxia/Hypoxia system designed for *in-vitro* cell culture experiments

Our Baker Grow solutions are packed with new, innovative features that allow you to study even the most complex cell interactions under perfect physiological oxygen conditions. Whether you're hoping to replicate the environment of blood vessels or lung tissue, Baker provides the best tools for the job.



Pho₂x Box included features

GAS CONTROLLER

- Controls and monitors O₂ and CO₂ levels with a single touchscreen
- Intuitive touchscreen with large font number display
- Small footprint (325mm wide x 298mm deep)
- Easy user set up, takes around 10 minutes
- Only needs N₂ and CO₂ cylinders for operation
- 2 year warranty for peace of mind

CULTURE CHAMBERS

Your choice of 4 types of cell culture chamber. Each is gas tight, has removable shelving and is easy to clean.

- Black (light reducing) - small and large
- Clear - small and large



TOUCHSCREEN CONTROL

- O₂ control (from 0.1% to 20.0% in 0.1% increments)
- CO₂ control (from 0.1% to 20.0% in 0.1% increments)
- Hypoxic cycling
- Temperature display (of Cell Culture Chamber)

AUDIBLE ALARMS

- Low gas (either CO₂ or N₂)

DATA LOG

- Up to 12 months data history
- One data set stored per minute, each data set comprises: Time, Date, O₂ (Set/Actual), CO₂ (Set/Actual)
- Stored on SD card provided




Figure 1. blot provided by PhD Daniela Mennerich, University of Oulu, Finland, showing proof of hypoxic conditions in the Pho₂x Box inside a standard CO₂ incubator.



CULTURE CHAMBER DIMENSIONS

		Small (Black or Clear) Chamber		Large (Black or Clear) Chamber	
		mm	inches	mm	inches
External dimensions	Width	335	14	355	14
	Height	173	6.8	238	9.4
	Depth	369	14.5	369	14.5
Internal dimensions	Width	335	13.2	335	13.2
	Height	140	5.5	205	8.1
	Depth	280	11	280	11
	Volume	13.1 Litres		19.2 Litres	
Workstation capacity-96 well plates	Number of plates (128mm x 86mm x 17mm)	8 plates of shelf, (4 stacks of 2 plates high on shelf)		24 plates, (4 stacks of 2 plates high per shelf, 3 shelves)	
Number of flasks (150mm x 80mm x 86mm)	Number of flasks (150mm x 80mm x 86mm)	6 (on shelf)		18 (6 on each shelf)	
PhO ₂ x Box System	Kgs/Lbs	5.5kg/12.1 lbs (chamber) & 5.5 kg/12.1 lbs (controller)		6.5 kg/14.3 lbs (chamber) & 5.5kg/12.1 lbs (controller)	

Usable Internal Volume: 13.1 litres (Small Chamber) | 19.2 litres (Large Chamber)

		mm	inches
Gas Controller External dimensions	Width	325	12.8
	Height	238	9.4
	Depth	298	11.7

PERFORMANCE DATA


	Gas Type		Coefficient of Variation (based on >700 readings)	
	O ₂ %	CO ₂ %	O ₂ %	CO ₂ %
Settings on Controller	0.1	5	<8%	<5%
	1	5	<5%	<5%
	5	5	<2%	<5%

PERFORMANCE DATA AT DIFFERENT GAS LEVEL SETTINGS

	0.1% O ₂ / 5.0% CO ₂	
	O ₂ %	CO ₂ %
Mean	0.1	5.0
%CV	6.7	2.2

	1% O ₂ / 5.0% CO ₂	
	O ₂ %	CO ₂ %
Mean	1.1	4.9
%CV	7.9	1.9

	5% O ₂ / 5.0% CO ₂	
	O ₂ %	CO ₂ %
Mean	5.1	4.9
%CV	1.3	1.6



Additional specifications

ACCESSORIES

- Anoxic/Anaerobic operation using palladium catalyst and special gas mix
- Single cable gland (must be factory fitted)
- Gas Sample Port
- Stand-alone O₂ meter for validation


POWER SPECIFICATIONS

- Voltage: 100 Volts to 240 Volts AC
- Power: 15 Watts
- 130kWh/Year (Based on 24 / 7 / 365)


When inside a workstation, the PhO₂x Box provides a second, independently controlled atmosphere for experimental flexibility.




Time to Set Point Large Culture Chamber O₂-5% CO₂-5%



30 Second Door Opening Recovery Large Culture Chamber O₂-5% CO₂-5%



5 Minute Door Opening Recovery Large Culture Chamber O₂-5% CO₂-5%





 Baker



bakerco.com

